

## Disposable [Plaster Bandage](#) for Arms and Legs Bone Fractures Fixed

### 1.Description

This [plaster bandage](#) is composed of gauze bandage starched by slurry and gypsum powder. It can be hardened in a short time after soaking by water. So it has strong ability to shape and good stability.

### 2.Features

- When the bandage is soaking by water, few gypsum loss.
- Harden time can be controlled: 2-5mins(ultrafast), 5-8mins(fast), 4-8mins(usual). It also can be made according to customer's request.
- strong hardness
- The drying time is fast, it can completely dry within 36 hours.
- Strong adaptability, high temperature resistant (+ 40 °C), cold (-40 °C) non-toxic, no stimulation, no irritation.

### 3.Application

Suitable for bone fractures fixed, deformity correction, inflammation limbs braking, osteomyelitis, bone tuberculosis, bone tumors and joint keratoplasty body fixed model and mould making, etc

### 4.Available Size & Package

Size	Package(rolls/ctn)	Meas(cm)
5×270cm	240	54×32×25
7.5×270cm	240	54×32×35.5
10×270cm	120	54×32×24
12.5×270cm	120	54×32×28
15×270cm	120	56×32×33
20×270cm	60	56×32×23

### 5.Notice

- Stored in the relative humidity is not more than 80%, no corrosive gas and well-ventilated room.
- Shelf life: 24 months



**Please feel free to contact us for some samples!**

**Contact: Lavender**

**Email:**[sales1@lanyuansupply.com](mailto:sales1@lanyuansupply.com)

**Tel&WhatsApp:**0086-13657299068