

[Factory Combat Emergency Tourniquet](#)

Description

- 1.Outdoor quick release tourniquet, one-hand operation.
- 2.Lightweight and portable, elastic and high strength.
- 3.Easy use, taut over perforation, quick hemostasis.
- 4.It can also be used for bundling large fixed object.
- 5.More suitable for outdoor survival, emergency medical us

Name	Tourniquet
Color	Black, Green,Blue,Red
Instrument classification	Class I
Material	Polyester+Aluminum+Nylon
Size	65cm, 75cm , 95cm, or Custom Size
Type	Stop Bleeding, First-Aid Devices
Design	Metal rod, with an iron clip
Application	First-aid
Function	Emerfergency Care
Place of Origin	China
Power Source	Manual
Feature	Portable
Logo	Accept Customized Logo
Certificate	CE , ISO
OEM	Avalibale

Features

1. Featuring NAR's Red Tip Technology™ with a red elliptical tip to assist user in locating and threading during application
2. Reinforced windlass clip and highly visible security tab includes a writeable area to record the time of application
3. Its effectiveness as a life-saving instrument for our soldiers on the battlefield has played a significant role in improvement of modern warfare combat injury medical care.
4. The Tourniquet has revolutionized battlefield care for the user, given its ISR (Institute of Surgical Research) confirmed effective function and its lightweight

Advantages

1. With anodized surface, the CNC machining aluminum handle has good maintenance of strength in extreme conditions. Black.

2. The two triangular rings is the dual locking system which keeps the fixation when movement is required.
3. The metal screw clamp, with tooth type clip and metal screw (anti-rust and with lubricant coating), provides more precise strength and fixation to avoid relaxing the tourniquet during movement.
4. Without self-adhering band design, Tactical Tourniquet can be used in dirty, muddy, sandy, also in the rain day and snow day.

Scenes To Be Used

Easy to Use: Rapid hemostatic tourniquet and easy to apply with one hand tourniquet on your on limbs

Life Saving: Awesome for first aid when hunting during camping and hiking,100% effective in blocking blood flow in both upper and lower limbs

Variable Usage: tourniquets are suitable for all situations

Easy to Carry: The compact size is easy to carry, it could be put in your backpack or pocket

Good Quality: Made of nylon and magic stick; Easy to clean and can be reused

[China emergency tourniquet supplier](#)

Use

- 1.Tactical Tourniquet is constructed of a solid piece of aircraft aluminum handle, two high impact plastic, acetyl triangular rings and a safety set screw.
- 2.With about 1m webbing band, it can be adjusted to fit any part of body or any age of persons.
- 3.It is chosen by many of the most elite and experienced warriors because of its outstanding ability to control severe bleeding and high level of reliability and ease of application.
- 4.2 or 3 pieces can be applied together as a bone fracture fixing band.

Payment and delivery

Package	1pc/bag, 100pcs/ctn
Payment	TT,L/C,Western Union
Delivery Time	Usually about 25 days after receipt of your payment or based on your order.
Nearest Port	Wuhan, Shanghai, Shenzhen









FAQ

Q1. Why Choose LanYuan Protective Co., Ltd?

1. Professional Manufacturer for 15 years with advanced facility
2. Competitive price
3. Accept Custom Order
4. Offer OEM Service
5. Fast service and on-time delivery

Q2: How to place sample?

Welcome to contact the sales to get the free samples.

Q3: What's your leading time?

2-5 working days for goods in stock after payment was done.

Q4: What shipping methods do you offer?

Express :DHL/FedEx/UPS/TNT/Aramex etc

Airfreight /Sea Freight/Rail freight or support shipment to your China agent.

Q5: Do you accept OEM/ODM?

Yes,we can offer OEM /ODM Service

Contact Us

H.K LanYuan Protective(Wuhan)Co.,Limited

Wuhan LanYuan Protective Co.,Ltd

Milan

M.& Whatsapp: 0086-18062501536

Tel: 0086-027-83734835

E-mail:milan@lanyuansupply.com

Please feel free to contact us for enquires and you will get reasonable prices and best quality.

[China tourniquet manufacturer](#)