

High Quality [Polyester Tubular Bandage](#) for Limbs and Trunks

Three types:

1. 100% polyester tubular bandage
2. 100% [Cotton tubular bandage](#)
3. High elastic tubular bandage

Benefits:

360 degree compression and support;
seamless, knit tubular construction;
ravel-free, run-free cutability;
secures without pins, clips or tape;
easy to cut for joint/appendage openings

Specifications:

Width	Length	Suggested Application
5cm	22.86m	Small limbs
7.5cm	22.86m	Medium/large limbs
10cm	22.86m	Extra large limbs
15cm	22.86m	Adult Trunks
20cm	22.86m	Large Trunks

Any width and length can be customized as your required.







Applications:

1. For use in treating tissue injuries, edema reductions, sports injuries, post-burn scarring, sprains, strains and for post-cast support, securing wound dressing and covering plaster bandages and dressings.
2. For use as a double layer for additional compression and support.

How to use?

Before use plaster of paris bandage or casting tape, according to size of injured part, choose width, cut suitable length, and set bandage on body part, then pull up.

Notice:

Remark: non sterile products, pls don't touch wound directly.

Storage: pls store in cool and dry area.

Racher

Whatsapp: 008615527106178

Skype&Email: sales2@lanyuansupply.com