Manufacturer Urine Collection Leg Bag

Product Description

Urine Bag is made from PVC in medical grade. It consists of bag, connecting tube, taper connector, bottom outlet and handle

It is intended to use with an indwelling catheter in persons who are incontinent of urine, can not urinate in the normal way, or need to have the bladder flow continually

Specification:

Product name:	Urine Collection Leg Bag
Bag Material:	Medical grade PVC
Bag Thickness:	1.2mm
Volume:	750ml
Tube Size:	Diameter:6.5mm Length:10cm or customized
Sterile:	sterile
Certificate:	CE&ISO13485
Attention :	for single use
Supply Ability:	1,000,000 pcs per day

Pictures Display: Medical supply disposable urine drainage bag

Advantages:

1. Professional Manufacturer for 10 years with advanced facility;

2. Easy to use and disposable

3.Factory price with best quality;

 $4.Accept \ Custom \ Order$, available in different sizes ;

5.Offer OEM Service;

6.Fast service and on-time delivery;

7. Main Market is Europe, America, Southeast Asia, etc.;

Packaging & Elivery:

Package:	10 pc/PE bag ,250 pcs/ctn
Deliver:	15-25 days or according the quantity









China medical leg bag supplier

FAQ

Q1:Whay Choose LanYuan Protective Co.,Ltd?
1.Professional Manufacturer for 15 years with advanced facility
2.Competetive price
3.Accept Custom Order
4.Offer OEM Service
5.Fast service and on-time delivery

Q2: How to place sample? Welcome to contact the sales to get the free samples. Q3: What's your leading time?2-5 working days for goods in stock after payment was done.

Q4: What shipping methods do you offer? Express :DHL/FedEx/UPS/TNT/Aramex etc Airfreight /Sea Freight/Rail freight or support shipment to your China agent.

Q5: Do you accept OEM/ODM? Yes,we can offer OEM /ODM Service

Contact Us

H.K LanYuan Protective(Wuhan)Co.,Limited Wuhan LanYuan Protective Co.,Ltd

Yuan (Sales department) M.& Whatsapp: 0086-18827431846 Tel: 0086-027-83734835 E-mail:sales1@lanyuansupply.com

We are the <u>China Urine Collection Leg Bag factory</u>, welcome to contact us for the samples !