

Orthopedic Cast Padding Wholesale at Factory Price

Three types:

1. 100% polyester tubular bandage
2. 100%cotton tubular bandage
3. High elastic [tubular bandage](#)

What the Applications of Medical Bandage:

1-for use in treating tissue injuries, edema reductions, sports injuries, post-burn scarring, sprains, strains and for post-cast support, securing wound dressing and covering plaster bandages and dressings.

2-for use as a double layer for additional compression and support.

What the Specification of our Tubular Bandage:

Width	Length	Suggested Application
5cm	22.86m or as Your Required	Small limbs
7.5cm	22.86m or as Your Required	Medium/large limbs
10cm	22.86m or as Your Required	Extra large limbs
15cm	22.86m or as Your Required	Adult Trunks
20cm	22.86m or as Your Required	Large Trunks

Any width and length can be customized as your required.







What the Benefits of our Surgical Bandage:

1. 360 degree compression and support;
2. Seamless, knit tubular construction;
3. Ravel-free, run-free cutability;
4. Secures without pins, clips or tape;
5. Easy to cut for joint/appendage openings

How to use?

Before use plaster of paris bandage or casting tape, according to size of injured part, choose width, cut suitable length, and set bandage on body part, then pull up.

Notice:

1. Remark: non sterile products, pls don't touch wound directly.
2. Storage: pls store in cool and dry area.

Payment and delivery:

Payment	TT,L/C,Western Union
Delivery time	Usually 25 working days after earnest or based on your quantity

Racher

Whatsapp: 008615527106178

Skype&Email: sales2@lanyuansupply.com