

Function:

1. Permanent spider vein removal
2. Capillaries expansion removal
3. Cherry hemangioma removal
4. Vascular therapy
5. All kinds of [Blood Vessel Removal](#)

Product parameters:

Item	Portable High-frequency Blood Vessel Removal Beauty Equipment
Operation System	High Frequency
Display	5.0 inch color intelligent touch screen
Treatment Needle	0.01mm diameter
Output Frequency	6.78MHZ
Therapy Handle Length	1.5M
Working Voltage	AC: 220V+10%, 50HZ; 110V+10%, 50HZ
Fuse Specification	Max 5A
Language	Chinese/English
Application	Beauty salon, clinic, hospital and home

Working Principles:

1. Vein Removal is a technologically advanced procedure designed to diminish leg and facial spider veins using a highly concentrated beam of RF energy. Spider vein removal and vein removal are popular because they are minimally invasive, virtually pain free, and require no general anesthesia.
2. In spider vein removal, High-RF passes through the skin and sends pulses of energy to the targeted vein without harming the surrounding tissue. As the blood absorbs the High-RF, it heats up and causes the targeted vein to shrink, eventually resulting in its disappearance. Spider vein removal requires no injections and no incisions. Generally, patients only experience a minimal amount of heat. Treatment time is about 20 minutes but may be longer or shorter depending on the severity of the veins and the size of the affected area. For optimal results, several spider vein removal sessions may be necessary.
3. With Spider vein removal, veins fade significantly or clear completely. Skin looks more uniform and youthful. Although Vein removal may not correct all the imperfections of the treated area, the overall results are successful; most patients are satisfied with treatment. Vein removal does not prevent new veins from forming, however. If new veins appear over time, further treatment may be needed to maintain the results of previous treatment.

Advantages:

1. It can continue to work for 24 hours.
2. Its output power is more stable, ensuring the treatment consistency and safety.
3. There are two working mode for doctors to choose: continuous mode and super pulse mode.
4. The professional handle ensure the protection of normal skin and the effect of couperose skin.

Packaging&delivery:

Packaging	1.Styrofoam 2.Fumigation wooden case 3.Air box
Delivery	Based on your order quantity

Product exhibition:

