

Supply Disposable Urinary Incontinence Bed Pee Pads

Products Description

Cottony topsheet, with rhombic embossment for quick absorbing , high quality fluff pulp and SAP mixed absorbent core ,bottom tissue paper layer for extra absorbency, the waterproof film backsheets guarantees hygiene and cleanness.

Name:	Disposable Underpads
Item NO.:	LY-UP-002
Material :	Non-woven fabric surface +Tissue paper +Fluff pulp +SAP+Tissue paper+PE film
Color :	Blue,Red, Green,White
Size:	40x60cm, 60x90cm, 80x180cm or Customized
Weight:	35~120g/pcs or as requested
Absorption:	200~3000ml or as order requested
Age Group:	all ages
Anti-Leak :	3D Leak Prevention Channel
Features:	Super Absorbent , Dry Surface ,Plain Woven ,Disposable, Eco-friendly, Comfortable
Place of Origin:	China

Material:

- * 1st Layer: Non Woven or Dry Mesh Cover
- * 2nd Layer: Tissue Paper
- * 3rd Layer: Fluff Pulp Mix Super Absorbent Polymer(SAP or GEL)
- * 4th Layer :Tissue Paper
- * 5th Layer: PE Film

Main characteristics:

1. Super absorption
2. Non-woven fabric
3. Soft and comfortable
4. Hygienic and clean

Advantage:

- 1.Upper layer of non-woven fabric for quick water absorption, refreshing and deodorizing.
- 2.It is disposable and healthy, one piece can be used for entire day.
- 3.Environmentally-friendly fabric, soft and breathable, tear-resistant and leak-proof for your

dog.

4.It will provide a great protection and cleanliness for your house.


WATERPROOF TEST



① Put a paper towel under the underpad



② Pour 200ml of water



③ Water is quickly absorbed



④ The bottom paper towel is still dry



⑤ Take a paper towel and press the surface directly



⑥ The surface of the paper towel is still dry

Features:

* **HEAVY ABSORBENCY** - Disposable underpad features a heavy absorbency core to provide incontinence protection for beds, furniture and floors

* **SOFT QUILTED TOP SHEET** - Disposable bed pad with a soft diamond quilted top sheet is

gentle on sensitive skin

* CORE ABSORBS FLUID AND LOCKS AWAY ODOR - The fluff and polymer core absorbs fluid and quickly locks away odor

* The large size of this underpad makes it great for use as a bed pad, furniture protector or puppy pee pad.

* CONVENIENTLY PACKAGED 50 PER CASE - These chux are pre-folded and can easily fit in a bag for on the go convenience.

[China Disposable Absorbent Underpads supplier](#)


Packaging Detail: 10pcs/bag, 20pcs/bag neutral package with plastic bags or colorful printing packing from customer's order and customer's logo

Suitable for: Hospital postpartum, postoperative, coma patient care and adult incontinence and infant care, women's menstrual care.

Scope: Applicable to hospital gynecology, obstetrics, orthopedic examination, operating table, B-, CT, endoscopy and other local health quarantine inspections.

Understand The Product

I care for you every day.







Package and Delivery:

Package :	10/15/20/25pcs/pack
Packing Type :	PE Bag /Colored Bag
Quantity./Container :	55000pieces / 20FT ;160000Piece /40HQ
Lead time :	10-25days from deposit
Payment Term :	TT , LC , West Union

FAQ:

Q1.Why Choose LanYuan Protective Co.,Ltd?

1.Professional Manufacturer for 15 years with advanced facility

- 2.Competitive price
- 3.Accept Custom Order
- 4.Offer OEM Service
- 5.Fast service and on-time delivery

Q2: How to place sample?

Welcome to contact the sales to get the free samples.

Q3: What's your leading time?

2-5 working days for goods in stock after payment was done.

Q4: What shipping methods do you offer?

Express :DHL/FedEx/UPS/TNT/Aramex etc

Airfreight /Sea Freight/Rail freight or support shipment to your China agent.

Q5: Do you accept OEM/ODM?

Yes,we can offer OEM /ODM Service

Contact us :

H.K LanYuan Protective(Wuhan)Co.,Limited

Wuhan LanYuan Protective Co.,Ltd

Yuan (Sales department)

M.& Whatsapp: 0086-18827431846

Tel: 0086-027-83734835

E-mail:sales1@lanyuansupply.com

If you have any questions , please always feel free to contact us

,we are the [China incontinence bed pee pads factory](#) , looking forward to your inquiry.