

# Lanyuan 808nm Diode Laser Hair Removal Machine

## Major function

### \*Hair Removal

Get rid of unwanted hair like lip hair, beard, armpit hair, hair on arms, legs and bikini area. Suitable for both light and dark skin.

### \*Skin Rejuvenation



## Working theory

1. The instrument can effectively penetrate to the targeted tissue, and the appropriate laser pulse time can ensure the accumulated energy to break the targeted tissue without damaging the surrounding tissues.
2. The appropriate energy density can provide enough energy output in a certain time to break the targeted tissues without damaging the surrounding tissues.
3. The appropriate epidermal protection measure can ensure enough damage of the targeted tissues without damaging the surrounding tissues;
4. The multi-pulse laser of this semiconductor laser hair removal series chooses the mode of a relative low energy density to heat mascara to 75°C, and maintain this state with the hand tools sliding under 10Hz, then both the hair follicle and stem cell growth will become inactive, reaching the goal of permanent hair removal .

|                               |                                       |
|-------------------------------|---------------------------------------|
| Item                          | 808nm Diode Laser Hair Removal        |
| Type                          | Laser                                 |
| Laser type                    | Semiconductor/freezing painless laser |
| Q-Switch                      | No                                    |
| Style                         | Stationary                            |
| Laser instrument working mode | Pulse mode                            |
| Cooling Fin                   | 240*2                                 |
| Wave length                   | 808nm±10nm                            |
| Light spot size               | 10mm*16mm                             |
| Pulse width                   | 10-300ms                              |
| Maximum energy of a pulse     | 3J-140J                               |
| Pulses energy                 | 140J+10J                              |
| Laser instrument power        | 600W                                  |
| Input power                   | 1200W                                 |
| Working voltage               | 110V+ 10%,50Hz                        |
| OSD Language                  | Chinese, English                      |

|                         |               |
|-------------------------|---------------|
| Insurance specification | Max 10 A      |
| Size (Machine)          | 970*730*355mm |
| Weight                  | ≈47kg         |

## Highlights

### Better customer experience

**Cool and painless:** Germany micro channel cooling technology realize freezing point hair removal, fully guaranteed the sapphire skin protection system keep at the temperature of 4°C.

**Quick and complete:** Large square spot increase the treatment area and covers all the corners of skin.

**No hyperpigmentation:** penetrate directly to the targeted tissue.

### Safer

**Skin temperature control:** the cooling technology ensures the skin in comfortable condition.

**No inflammation:** the laser energy reach directly to the hair follicles without damaging the skin.

**Latest technology:** guaranteed by the latest technology in Optics Valley in China, where lies the high-density high-tech enterprises.

### More efficient

**Large spot:** 12mm\*16mm big spot make it more easier for large area hair removal.

**User friendly:** multi-language menu suitable for different operators

**Accurate:** High energy density output ensure it suit for all kinds of skin. Laser energy can reach deep hair follicles and be effectively absorbed.

**Flexible:** adjustable frequency, maximum output frequency is 10HZ

### Longer service life

**Cooling system:** Air + Water + Diode + TEC Cooling. TEC is the latest cooling method which widely use in refrigerator, this new cooling method can confirm the diode laser in a more suitable working environment and control it in low temperature even for long time contentiously working.

**Strict temperature control:** if the temperature is over 35 degree, the machine will stop work, this can well protect the therapy handle safety, especially for the bars and crystal filter.

**Imported power supply:** Taiwan imported, the best power supply in the market, make sure the machine stable and energy enough.

**Germany imported bars:** confirm the shoting of the machine over 5,000,000 times and have enough energy.

### Lower running cost

With less treatment time and longer service life, the unit running cost on each customer is considerably reduced.

## In Lanyuan, you can find the most suitable machine for your beauty career.

Customized function combination to meet different beauty purpose.

Various cases options to match different style in your work room.

### Service:

1. Technical instruction documents.
2. Technical support by pictures, videos and online instruction.
3. Response within 24 hours.

### Warranty:

One year warranty, life-long maintenance. Our engineers' technical support will be constant.

## Package

All instruments well packed by styrofoam, air box or fumigation wooden cases.



## Shipping and delivery

Shipping: by air, by sea or by express.

Delivery Time: 3-10 business days by sample.

**H.K LanYuan Protective(Wuhan)Co.,Limited**  
**Wuhan LanYuan Protective Co.,Ltd**

Flora Xie

Tel: 0086-027-83734835

Skype: jyboard

E-mail: [flora@lanyuansupply.com](mailto:flora@lanyuansupply.com)