

Functions:

1. [Wrinkle removal](#)
2. Skin rejuvenation
3. Acne treatment
4. Pigment removal
5. Blood vessel removal
6. Hair removal
7. Skin firming and reform, black doll skin rejuvenation
8. Slimming
9. Striae gravidarum improving
10. Freckle removal
11. Tattoo removal, nevus removal
12. Whiten and beautify skin

Product descriptions:

Type	Multi-Function Beauty Equipment
Treatment heads	4 heads
Spot area	10*50mm
Interval of light	2~5s
input consumption	≤1200w
Power	800w
Capacitance	30000UF
Laser type	Passive Q Switch ND YAG laser
Laser Bar	Diameter 5mm*85mm
Laser Driving power	500w
Energy density	1~50mj/cm ²
Total pulse width	150ms
Pulse interval	20ms
Cooling Fin	Single Fin
Operation panel	10.4 inch
Size (Machine)	970*730*355mm
OSD Language	Chinese, English, and Korean
Material	ABS + Aluminium Alloy

Advantages:

1. Combine E- light and ND YAG Laser.
2. Separate Power Supply, E-light and ND YAG Laser don't effect each others.
3. E-light use the latest OPT technology, UK Firstlight Lamp with at least 300,000 times and Artificial Sapphire filter crystal.
4. ND YAG Laser with three working probes, 1064nm and 532nm for pigment removal, 1320nm balck doll carbon treatment for skin rejuvenation.

Packaging&delivery:

Packaging	1.Styrofoam 2.Fumigation wooden case 3.Air box
Delivery	Based on your order quantity

Product exhibition:







Laser with hand



The laser head



Elight working Handle



RF working handle



RF working head