

Therapy scope:

- 1,Wrinkle removal,acne removal,freckles,freckle spot
- 2,Dispel black rim,hemangioma, nevus flammeus
- 3.[Skin rejuvenation](#),skin improvement
- 4.Pore shrink,telangiectasia
- 5.Hair removal

Technical parameters:

Item	OPT-E Light Perfect Pulse Beauty Equipment
Size	1180*500*385mm
Treatment Handle	Single handle tool(532nm-950nm or 640nm- 950nm)
Light spot area	10mm*50mm
Energy density	1-50J/cm2,stepping 1J Light
interval	2-5s(adjustable),stepping 1s Input
power	Max 1200W
Material	ABS + Aluminium Alloy
RF power	Max 100W
Cooling method	Circulation cooling with the internal distilled water,air cooling and semiconductor refrigeration
Language	Chinese,,English and Korean

Packaging&delivery:

Packaging	styrofoam, air box or fumigation wooden cases
Delivery	Based on your order

Details:



523nm handle tool



640nm handle tool

Factory:

