

Functions:

1. [Scar removal](#), especially at hypertrophic scar, pock, and striae gravidarum
2. Stretch mark removal
3. Wart and mole removal
4. Hemangioma removal
5. [Wrinkle removal](#)
6. Skin rejuvenation, make skin fine and smooth, resilient

Product descriptions:

Type	Laser
Laser Type	CO2 Laser
Q-Switch	No
Style	Stationary
Light guide system	Seven joint arm
Laser power	30W
Wavelength	10600nm
Aiming wavelength	635nm
Pulse width	0.1ms-10ms(adjustable)
Pulse interval	1ms-5000ms (adjustable)
Pulses energy	0.1-300mj (adjustable)
Focal spot diameter	0.1mm--2mm (adjustable)
Input power	600W
Cooling system	Air cooling, water cooling
Operation system	10.4 inch color intelligent touch screen
OSD Language	Chinese, English
Chassis material	ABS environmental friendly medical chassis

Advantages:

1. Combine the conventional therapy of pulse laser with the fractional laser technology, and these two modes can be shifted in 2 seconds.
2. The treatment is painless and its time is short, besides, there is no need for additional consumables.
3. By promoting the collagen hyperplasia and reorganization, its treatment effect can be kept for several years.
4. The minimum diameter of the laser is 0.075mm, so its damage is minor and one can recover quickly.
5. As the main treatment area remains intact during operation, one can use cosmetics at the second day after the operation.

Packaging&delivery:

Packaging	1. Styrofoam 2. Fumigation wooden case 3. Air box
-----------	---

Delivery	Based on your order quantity
----------	------------------------------

Product exhibition:





Contact: Racher Cheng

Tel.& Whatsapp: 0086-15527106178

Skype: sales2@lanyuansupply.com

Email: sales2@lanyuansupply.com