Functions:

- 1. Vaginal Tightening
- 2. Vaginal Beauty
- 3. Vaginal Lubrication
- 4. Vaginal Protecting

Product datas:

Item Vaginal Tightening Fractional Co2 Laser Salon

<u>Equipment</u>

Type Laser

Laser Type CO2 Laser

Q-Switch No

Light guide system Seven joints(Korea)

Light power 30w(40w)

Output Single point, Continuous, super pulse, dot matrix

Cutting head(50mm,100mm) 2+4

Wavelength 10600nm

Vagina Tightening Laser protective jacket Titanium alloy protective jacket*2

Cooling Fin Double fin

Dot distance 0.1-2.6mm (adjustable)
Pulses energy 0.1-180mj (adjustable)

Focal spot diameter 100um-2000um

Focal spot quantity 1(1*1)-40000(200*200) (adjustable)

Figure square 0.1mm*0.1mm-20mm*20mm(adjustable)
Working voltage AC:220V+ 10%, 50Hz; 110V+ 10%, 50Hz

Screen 10.4 Inch Color Touch Screen

Color Silver
Weight ≈42kg

Advantages:

- 1. Minimally invasive, virtually painless procedure and fast recover, private laser is painless, simple, and safe
- 2.It can cure urinary incontinence, vaginal relaxation, the darkness of the private part, and i nflammation
- 3.Its effect can maintain1-1.5 year with a course of treatment
- 4.It has no trauma and pains, even without anesthesia
- 5. After operation, it leaves no discomfort, and it needs no special care
- 6. Allow having sexual life within 10 days
- 7.It only needs 20 minutes to complete a course of treatment
- 8.It has no side effects with repeated therapy

Packaging&delivery:

Packaging	export standard wooden box and polyfoam
_{ll} i uckaging	caport standard wooden box and poryroun

Product exhibition:







